

## SPACE DOSIMETRY WITH A 3D SILICON DETECTOR TELESCOPE – THE ISS VERSIONS OF TRITEL

Attila Hirn, Tamás Pázmándi, Sándor Deme, István Apáthy, Antal Csőke, László Bodnár\*

Hungarian Academy of Sciences KFKI Atomic Energy Research Institute \*BL-Electronics

hirn@sunserv.kfki.hu

10-12 September, 2007

12th WRMISS

#### Outline



- ISS versions of TriTel
  - 3D silicon detector telescope
  - Measurement ranges of the instrument
  - Construction
- Electron sensitivity of the Pille dosimeters (off-topic)
- Summary

# The TriTel 3D telescope



Z telescope

X telescope

10-12 September, 2007

12th WRMISS

- 3 x 2 db 300 µm Si detector
- Detector chips of the telescope are connected as AND gate in coincidence
- D /from dep. energy spectra/
- LET spectra and Q /from the gated spectra/

→ dose equivalent (H)

- 3D sensitivity
- r = 6.9 mm
- p = 8.7 mm



## LET and $\Delta E$ ranges

- Minimum deposited energy in the coincidence spectra:
  - By relativistic protons (LET<sub>min,Si</sub> = 0.4 keV/ $\mu$ m), 0° incidence
  - $-\Delta E_{\min, \text{ classical}} = LET_{\min} * w = 120 \text{ keV}$
- But thin detectors!  $\rightarrow$  Landau distribution
- Most probable energy loss:



$$\Delta_0 = \xi \left[ \ln \frac{3000\beta^2 \xi}{Z^2 (1 - \beta^2)} \right] (eV); \text{ where } \xi = \frac{1.54 \cdot 10^5 \rho}{\beta^2} \cdot \frac{Z}{A}$$
  
and  $FWHM = \frac{4.0\xi}{\Delta_0}$ 

12th WRMISS

### LET and $\Delta E$ ranges



Minimum @ 2-5 GeV

$$\Delta E_{\min} = \sim 50-60 \text{ keV}$$

- $LET_{max,Si} = 230 \text{ keV/}\mu\text{m}$ 
  - (120 keV/µm in water)
- $\rightarrow \Delta E_{\text{max}} = \sim 85 \text{ MeV}$
- $\Delta E > 85 \text{ MeV} \rightarrow$  overflow bin

# TriTel measurement ranges

Physical quantity	Measurement (signal) ranges
LET in Si	0.4 keV/μm – 230 keV/μm
Energy imparted to the detector material	60-80 keV – 85 MeV
Charge created by an incoming particle	2.7-3.6 fC – 3.8 pC
Voltage step	1.3-1.8 mV – 1.9 V
Amplitude of the signal to be processed by the ADC	2 mV – 2 V
	Physical quantity   LET in Si   Energy imparted to the detector material   Charge created by an incoming particle   Voltage step   Amplitude of the signal to be processed by the ADC



AEKI

- 12-bit ADC
- LUT  $\rightarrow$  4096 channel spectra
- 3 x 2 x 2 deposited energy spectra per day
- 3 time spectra
   (10 s resolution
   + SAA flag)





# Construction of the detector unit



### User interface unit



- PC104 module with graphical display, 6 menu buttons, 1 Ethernet and 2 USB ports
- Parameters can be changed manually
- Data written to memory card
- Displaying the results of the preliminary data analysis
- Possible software updates

## **SURE TriTel**



- Upload expected in 2009 (Increment 19)
- Operation for 6 months
- Comprises
  - TriTel 3D silicon detector telescope
  - Palfalvi's passive detector stacks (3 stacks with axes coincident with the axes of TriTel) consisting
    - several layers of SSNTDs (plastics and converter layers) and
    - Al<sub>2</sub>O<sub>3</sub> TL detectors (integral absorbed dose measurement)

10-12 September, 2007





 $25 \ \text{mm} \ \text{x} \ 55 \ \text{mm} \ \text{x} \ 20 \ \text{mm}$ 

## **ISS TriTel**



- In the Russian segment of the ISS
- In cooperation with IBMP
- Upload expected in 2009
- Operation for at least 6 months
- Comparing the results with SURE TriTel
- Later on  $\rightarrow$  external unit  $\rightarrow$





## Electron sensitivity of Pille during EVAs (1)





- EVA → astronauts are protected by the Orlan-M spacesuit
- $\sim 1.2 \text{ mm Al(eq.)}$  shielding
- Pille dosimeter:
  - Vacuum bulb (glass)  $\sim 1.09 \text{ mm Al}$
  - Cylindrical Al holder  $\sim 0.75$  mm Al
  - Stainless steel tube  $\sim 0.71 \text{ mm Al}$
  - Protective Al case ~ 0.75 mm Al
  - Polyamid insert ~ 0.89 mm Al
  - $\sim$  4.2 mm Al(eq.) shielding provided by the different parts of the Pille dosimeter
- $\rightarrow$  considerable fraction of the electrons are shielded  $\rightarrow$  underestimate of the electron dose

## Electron sensitivity of Pille during EVAs (2)

- Substituting the protective case (and the polyamid insert) with a thin Kapton foil (~ 0.05 mm Al)
- $\rightarrow$  2.6 mm Al(eq.) shielding
- Expected doses behind the shielding were estimated using the Geant4 based Multi-Layered Shielding Simulation Software (MULASSIS)
- Simulations were run in planar geometry
- Proton and electron spectra taken from AP8 and AE8 model for ~ISS orbit (altitude of 400 km and an inclination of 51.6°) [Redell and Atwell, 2005]

## Electron sensitivity of Pille during EVAs (3)



Redell & Atwell

## Electron sensitivity of Pille during EVAs (4)

• Absorbed doses in skin tissue relative to doses behind Orlan-M:

	D/D <sub>Orlan</sub> (0°, parallel)	D/D <sub>Orlan</sub> (cosine law)
Orlan-M	100%	100%
Pille with metal case	5%	4%
Pille without metal case	21%	18%

# Electron sensitivity of Pille during EVAs (5)

• Absorbed doses in skin tissue relative to doses behind Orlan-M (cont'd):

0°parallel	D <sub>electr</sub> /D <sub>Orlan</sub>	<b>D</b> <sub>proton</sub> / <b>D</b> <sub>Orlan</sub>	D <sub>total</sub> /D <sub>Orlan</sub>
Orlan-M	4.5%	95.5%	100%
Pille w/o metal case	0.9%	66.6%	67.6%
Pille w/ metal case	0.2%	55.0%	55.3%

Cos law	<b>D</b> <sub>electr</sub> / <b>D</b> <sub>Orlan</sub>	<b>D</b> <sub>proton</sub> / <b>D</b> <sub>Orlan</sub>	<b>D</b> <sub>total</sub> / <b>D</b> <sub>Orlan</sub>
Orlan-M	2.4%	97.6%	100%
Pille w/o metal case	0.4%	72.5%	72.9%
Pille w/ metal case	0.1%	58.7%	58.8%

#### Summary



- Measurement range of TriTel in terms of deposited energy will be 60 keV – 85 MeV
- $\Delta E > 85 \text{ MeV} \rightarrow \text{overflow bin}$
- SURE TriTel: TriTel + 3 passive detector stacks
  - LET spectra (0.2 keV/µm 120 keV/µm in water), absorbed dose, average Q, dose equivalent (separately in SAA)
  - Integrated absorbed dose for the mission (TL)
  - LET spectra, and directional information (SSNTDs)
- ISS TriTel in the Russian Segment of ISS
- Pille electron sensitivity can be increased by a factor of ~ 4 5 by removing the protective metal case

## Acknowledgements

- The SURE project "New Complex Method for Determining the Equivalent Dose of Astronauts" (contract number: RITA 026069) is funded by the European Commission within the sixth framework programme for RTD activities (2002-2006), Specific Programme "Structuring ERA" support to Research Infrastructures and within the framework of the C98057 PECS agreement. The authors would like to acknowledge the services contributed by the representatives of ESA, hereby.
- The project ISS TriTel is carried out in cooperation with the IBMP and with the financial support of the Hungarian Space Office.



### Thank you for your attention

hirn@sunserv.kfki.hu